

If you encounter errors, call the service desk and provide them with this information:

- Identify yourself and advise them that this is for a **Telestroke Urgent Clinical Consultation**.
- Provide your location and the VCU # (found on the sticker on the monitor and/or equipment).

Home screen

Tap the screen to wake up the system

Tap to open the **Control panel** (camera, selfview, etc...)



Tap to make a **Call**.

Volume controls
(Mute button is at bottom of screen during call).

⚠ The Home screen should always be displayed prior to connecting to a site; if not, this may mean that a previous call was not disconnected. If you do not see the Home screen, press **End** to disconnect prior to connecting to your site.

Camera control & selfview


Before a call

Tap **Control panel** to invoke the camera adjustment menu.




During a call

Tap **Selfview** or **Camera**.




Turn **Selfview** On/Off or make it fullscreen.

Change **Selfview** position.



Tap **Camera control** to adjust your camera using the zoom +/- and direction arrows.



Receiving the call

- The Neurologist will connect to you and you will see them on the monitor once they have connected.
- Press the **Microphone** key to unmute and ensure other site can hear you.
- Assist with the patient assessment throughout the session by moving the camera or the patient as requested by the Neurologist.
- Press the red **End Call** button on the remote to disconnect at the end of every call.